



# BI Sobriotor®

A multi-language alcohol monitoring device with biometric voice verification technology.

## Key Features

### Compatibility

- BI HomeGuard® 200
- BI ExacuTrack®
- BI VoiceID®

### Battery Back-up

- Up to 12 hours
- Up to 4 tests

### Fuel Cell

- Calibration required every 6 months

### Language

- English
- Spanish

### Testing Schedule

- Regular intervals
- Random intervals
- On-Demand

## Remote Alcohol Monitoring Made Simple

Every year there are more than 1.5 million DUI arrests and approximately 3 million alcohol-related crimes committed. With the BI Sobriotor®, an offender who commits alcohol related offenses can be sentenced to monitored sobriety.

### Product Overview

BI Sobriotor is a hand-held, remote alcohol monitoring device that detects the presence of alcohol through an individual's deep lung breath sample. Sobriotor incorporates biometric voice verification technology, fuel cell alcohol testing technology, and BI's proprietary security technology to ensure the client is accurately identified and tested, the unit's integrity remains intact, and results are reported promptly to the host monitoring computer.

Monitoring with Sobriotor is performed remotely using a telephone connection, so that compliance with alcohol restrictions can be verified daily. The system allows officers to customize an alcohol breath testing schedule that fits the client's schedule. Test results are reported to the host monitoring computer and officers are notified according to agency notification procedures for results outside specified parameters.

### Easy Enrollment

The Officer Control Unit makes offender enrollment simple and fast. The unit's visual display and keypad guide the officer through client enrollment, enabling the officer to verify voice enrollment and perform on-demand field sobriety tests.

### Simple, Reliable Testing

To make remote testing consistently reliable, voice prompts are incorporated into Sobriotor that provide the offender with clear instructions. For example, if the individual is speaking too softly, Sobriotor prompts the offender to "speak louder." Voice prompts ensure that the offender is instructed through each step of the test.

### Officer-controlled Security

Variable security levels give the officer the ability to adjust the voice verification threshold for each offender. This feature reduces false alerts and minimizes attempts to circumvent the system. In addition, Sobriotor records the last failed voice verification test, which officers can access for voice comparison purposes.

### Sobriotor and Radio Frequency Monitoring: A Powerful Combination

Sobriotor may be used as a stand-alone sanction or in conjunction with radio frequency monitoring. In stand-alone mode, Sobriotor tests the offender randomly throughout the test period, as specified by the officer. When used with the BI HomeGuard® 200, the offender may be tested immediately upon entering the house, and randomly thereafter. In either case, only one phone line and one power outlet are required.

Monitoring of Sobriotor tests is performed at BI Monitoring Operations. Trained BI personnel perform testing procedures specified by the agency, and report all results. Agencies may also use an on-site host system allowing them internal control for compliance monitoring. BI works closely with every customer to determine which monitoring system best meets their needs.



**Call 800.701.5171 today to schedule a product demo.**